

C1R 抗原（重组蛋白）

中文名称： C1R 抗原（重组蛋白）

英文名称： C1R Antigen (Recombinant Protein)

别名： complement C1r subcomponent

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 506-705 amino acids of human C1R

技术规格

Full name:	complement C1r subcomponent
Swissprot:	P00736
Gene Accession:	BC035220
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	C1r B chain is a serine protease that combines with C1q and C1s to form C1, the first component of the classical pathway of the complement system.

mlbio 梅联生物
Good elisakit producers