

C2CD2 抗原（重组蛋白）

中文名称：C2CD2 抗原（重组蛋白）

英文名称：C2CD2 Antigen (Recombinant Protein)

别名：C2 calcium-dependent domain containing 2; TMEM24L; C21orf25; C21orf258

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 189-388 amino acids of human C2CD2

技术规格

Full name:	C2 calcium-dependent domain containing 2
Synonyms:	TMEM24L; C21orf25; C21orf258
Swissprot:	Q9Y426
Gene Accession:	BC062323
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	C2CD2 is an 696 amino acid protien that is secreted. The C2CD2 gene product has been provisionally designated C2CD2 pending f urther characterization.