

C12orf10 抗原(重组蛋白)

中文名称: C12orf10 抗原(重组蛋白)

英文名称: C12orf10 Antigen (Recombinant Protein)

别 名: chromosome 12 open reading frame 10; MYG; MYG1; Gamm1; MST024; MSTP024

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 156-376 amino acids of human C12orf10

技术规格

Full name:	chromosome 12 open reading frame 10
Synonyms:	MYG; MYG1; Gamm1; MST024; MSTP024
Swissprot:	Q9НВ07
Gene Accession:	BC028904
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	C12orf10, also known as MYG1, MYG1 is a 376 amino acid nucleo-mitochondrial protein belonging to the UPF0160 (MYG1) family. M YG1 is encoded by a gene that maps to human chromosome 12q 13.13 and is ubiquitously expressed in simple as well as complex e ukaryotes, with highest levels in testis. Considered to have a metal -dependent protein hydrolase (UPF0160) domain, MYG1 exhibits a mitochondrial targeting signal in the N-terminal region and a Pat7 -type nuclear localization signal in the region between amino acids



33-39. Although MYG1 displays differential patterns and levels of e xpression during embryonic development, expression in normal ad ult tissues is stable, suggesting MYG1 involvement in early develop mental processes and in adult stress/illness conditions. Elevation of MYG1 expression may be also associated with vitiligo susceptibility.