

C15orf40 抗原(重组蛋白)

- 中文名称: C15orf40 抗原 (重组蛋白)
- 英文名称: C15orf40 Antigen (Recombinant Protein)
- 别 名: chromosome 15 open reading frame 40
- 相关类别: 抗原
- 储 存: 冷冻 (-20℃)

概述

技术规格

Full name:	chromosome 15 open reading frame 40
Swissprot:	Q8WUR7
Gene Accession:	BC019820
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Encoding more than 700 genes, chromosome 15 is made up of approximately 106 million base pairs and is about 3% of the hu man genome. Angelman and Prader-Willi syndromes are associat ed with loss of function or deletion of genes in the 15q11-q13 region. In the case of Angelman syndrome, this loss is due to i nactivity of the maternal 15q11-q13 encoded UBE3A gene in the brain by either chromosomal deletion or mutation. In cases of P rader-Willi syndrome, there is a partial or complete deletion of t his region from the paternal copy of chromosome 15. Tay-Sachs disease is a lethal disorder associated with mutations of the HE



XA gene, which is encoded by chromosome 15. Marfan syndrom
e is associated with chromosome 15 through the FBN1 gene. Th
e C15orf40 gene product has been provisionally designated C15
orf40 pending further characterization.