

CA6 抗原（重组蛋白）

中文名称： CA6 抗原（重组蛋白）

英文名称： CA6 Antigen (Recombinant Protein)

别名： carbonic anhydrase 6; CA-VI; GUSTIN

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 18-179 amino acids of human CA6

技术规格

Full name:	carbonic anhydrase 6
Synonyms:	CA-VI; GUSTIN
Swissprot:	P23280
Gene Accession:	BC034350
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is one of several isozymes of carbonic anhydrase. This protein is found only in salivary glands and saliva and protein may play a role in the reversible hydration of carbon dioxide though its function in saliva is unknown.