CASK 抗原(重组蛋白)

中文名称: CASK 抗原(重组蛋白)

英文名称: CASK Antigen (Recombinant Protein)

别 名: calcium/calmodulin dependent serine protein kinase; CMG; FGS4; LIN2; TNRC8; CAGH39; CAMGUK; MICPCH; MRXSNA

储存: 冷冻(-20℃)

相关类别 抗原

概述

Fusion protein corresponding to a region derived from 727-926 amino acids of human CASK

技术规格

| Full name: | calcium/calmodulin dependent serine protein kinase |
|--------------------|---|
| Synonyms: | CMG; FGS4; LIN2; TNRC8; CAGH39; CAMGUK; MICPCH; MRXSNA |
| Swissprot: | O14936 |
| Gene Accession: | BC117311 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | This gene encodes a calcium/calmodulin-dependent serine protein ki nase. The encoded protein is a MAGUK (membrane-associated guan ylate kinase) protein family member. These proteins are scaffold pro teins and the encoded protein is located at synapses in the brain. Mutations in this gene are associated with FG syndrome 4, mental r etardation and microcephaly with pontine and cerebellar hypoplasia, and a form of X-linked mental retardation. Multiple transcript varian ts encoding different isoforms have been found for this gene. |