

小鼠抗 TSC22D4 单克隆抗体

中文名称：小鼠抗 TSC22D4 单克隆抗体

英文名称：Anti-TSC22D4 mouse monoclonal antibody

别名：TSC22 domain family member 4; THG1; THG-1; TILZ2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：TSC22D4

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	TSC22D4 is a member of the TSC22 domain family of leucine zipper transcriptional regulators (see TSC22D3; MIM 300506) (Kester et al., 1999 [PubMed 10488076]; Fiorenza et al., 2001 [PubMed 11707329]).[supplied by OMIM, Mar 2008]
Applications:	WB
Name of antibody:	TSC22D4
Immunogen:	Fusion protein of human TSC22D4
Full name:	TSC22 domain family member 4
Synonyms:	THG1; THG-1; TILZ2
SwissProt:	Q9Y3Q8
WB Predicted band size:	41 KD
WB Positive control:	Hela, K562 cell lysates



WB Recommended dilution:	500-2000
---------------------------------	----------