

CCNK 抗原（重组蛋白）

中文名称：CCNK 抗原（重组蛋白）

英文名称：CCNK Antigen (Recombinant Protein)

别名：cyclin K; CPR4

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 155-354 amino acids of human CCNK

技术规格

Full name:	cyclin K
Synonyms:	CPR4
Swissprot:	O75909
Gene Accession:	BC015935
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the transcription cyclin family. These cyclins may regulate transcription through their association with and activation of cyclin-dependent kinases (CDK) that phosphorylate the C-terminal domain (CTD) of the large subunit of RNA polymerase II. This gene product may play a dual role

le in regulating CDK and RNA polymerase II activities.