

CD44 抗原(重组蛋白)

中文名称: CD44 抗原(重组蛋白)

英文名称: CD44 Antigen (Recombinant Protein)

别 名: CD44 molecule (Indian blood group); IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; HCELL; HUTCH-I; ECMR-III

储存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 225-500 amino acids of human CD44

技术规格

Full name:	CD44 molecule (Indian blood group)
Synonyms:	IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; HCELL; HU TCH-I; ECMR-III
Swissprot:	P16070
Gene Accession:	BC004372
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system :	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a cell-surface glycoprotein involve d in cell-cell interactions, cell adhesion and migration. It is a receptor f or hyaluronic acid (HA) and can also interact with other ligands, such as osteopontin, collagens, and matrix metalloproteinases (MMPs). This



protein participates in a wide variety of cellular functions including lym
phocyte activation, recirculation and homing, hematopoiesis, and tumor
metastasis. Transcripts for this gene undergo complex alternative splici
ng that results in many functionally distinct isoforms, however, the full
length nature of some of these variants has not been determined.