

CDC45 抗原(重组蛋白)

- 中文名称: CDC45 抗原(重组蛋白)
- 英文名称: CDC45 Antigen (Recombinant Protein)
- 别名: cell division cycle 45; CDC45L; CDC45L2; PORC-PI-1
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 367-566 amino acids of human CDC45

技术规格

Full name:	cell division cycle 45
Synonyms:	CDC45L; CDC45L2; PORC-PI-1
Swissprot:	O75419
Gene Accession:	BC006232
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene was identified by its strong similar ity with Saccharomyces cerevisiae Cdc45, an essential protein require d to the initiation of DNA replication. Cdc45 is a member of the hig hly conserved multiprotein complex including Cdc6/Cdc18, the minich romosome maintenance proteins (MCMs) and DNA polymerase, whic h is important for early steps of DNA replication in eukaryotes. This protein has been shown to interact with MCM7 and DNA polymeras



e alpha. Studies of the similar gene in Xenopus suggested that this protein play a pivotal role in the loading of DNA polymerase alpha onto chromatin. Alternate splicing results in multiple transcript varian ts.