

CDHR4 抗原（重组蛋白）

中文名称： CDHR4 抗原（重组蛋白）

英文名称： CDHR4 Antigen (Recombinant Protein)

别名： cadherin-related family member 4; CDH29; PRO34300

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 18-180 amino acids of human CDHR4

技术规格

Full name:	cadherin-related family member 4
Synonyms:	CDH29; PRO34300
Swissprot:	A6H8M9
Gene Accession:	BC132756
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types (By similarity).

mlbio 梅联生物
Good elisakit producers