

CDKN2A 抗原

中文名称: CDKN2A 抗原

英文名称: CDKN2A Antigen

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Synthetic peptide corresponding to a region derived from 144-156 amino acids of human cyclin-dependent kinase inhibitor 2A

技术规格

Background	This gene generates several transcript variants which differ in their f irst exons. At least three alternatively spliced variants encoding distinct proteins have been reported, two of which encode structurally related isoforms known to function as inhibitors of CDK4 kinase. The remaining transcript includes an alternate first exon located 20 Kb upstream of the remainder of the gene; this transcript contains an alternate open reading frame (ARF) that specifies a protein which is structurally unrelated to the products of the other variants. This ARF product functions as a stabilizer of the tumor suppressor protein p5 as it can interact with, and sequester, the E3 ubiquitin-protein lig ase MDM2, a protein responsible for the degradation of p53.
Gene Accession	NP_000068
Swissprot	P42771
Synonyms	ARF; MLM; P14; P16; P19; CMM2; INK4; MTS1; TP16; CDK4I; CDKN2; INK4A; MTS-1; P14ARF; P19ARF; P16INK4; P16INK4A; P16-INK4A



		2
Full name	cyclin-dependent kinase inhibitor 2A	