

## CFAP298 抗原（重组蛋白）

中文名称： CFAP298 抗原（重组蛋白）

英文名称： CFAP298 Antigen (Recombinant Protein)

别名： cilia and flagella associated protein 298; Kur; FBB18; CILD26; C21orf48; C21orf59

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 1-290 amino acids of human CFAP298

### 技术规格

<b>Full name:</b>	cilia and flagella associated protein 298
<b>Synonyms:</b>	Kur; FBB18; CILD26; C21orf48; C21orf59
<b>Swissprot:</b>	P57076
<b>Gene Accession:</b>	BC000709
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that plays a critical role in dynein arm assembly and motile cilia function. Mutations in this gene result in primary ciliary dyskinesia. Naturally occurring readthrough transcription occurs from this locus to the downstream t-complex 10 like (TCP10L) gene. [provided by RefSeq, Apr 2017]

**mlbio** 梅联生物  
Good elisakit producers