

## CFB 抗原(重组蛋白)

中文名称: CFB 抗原(重组蛋白)

英文名称: CFB Antigen (Recombinant Protein)

别 名: complement factor B; BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf;

ARMD14

相关类别: 抗原

储 存: 冷冻(-20℃)

## 概述

Fusion protein corresponding to a region derived from 615-764 amino acids of human CFB

## 技术规格

Full name:	complement factor B
Synonyms:	BF; FB; BFD; GBG; CFAB; CFBD; PBF2; AHUS4; FBI12; H2-Bf; ARMD14
Swissprot:	P00751
Gene Accession:	BC004143
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes complement factor B, a component of the alter native pathway of complement activation. Factor B circulates in the blood as a single chain polypeptide. Upon activation of the alternat ive pathway, it is cleaved by complement factor D yielding the non catalytic chain Ba and the catalytic subunit Bb. The active subunit B b is a serine protease which associates with C3b to form the altern



ative pathway C3 convertase. Bb is involved in the proliferation of preactivated B lymphocytes, while Ba inhibits their proliferation. Thi s gene localizes to the major histocompatibility complex (MHC) clas s III region on chromosome 6. This cluster includes several genes i nvolved in regulation of the immune reaction. Polymorphisms in thi s gene are associated with a reduced risk of age-related macular d egeneration. The polyadenylation site of this gene is 421 bp from t he 5' end of the gene for complement component 2.