

## CLIC1 抗原（重组蛋白）

中文名称： CLIC1 抗原（重组蛋白）

英文名称： CLIC1 Antigen (Recombinant Protein)

别名： chloride intracellular channel 1; G6; NCC27

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 47-241 amino acids of human CLIC1

### 技术规格

<b>Full name:</b>	chloride intracellular channel 1
<b>Synonyms:</b>	G6; NCC27
<b>Swissprot:</b>	O00299
<b>Gene Accession:</b>	BC064527
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Chloride channels are a diverse group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, trans epithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 1 is a member of the p64 family; the protein localizes principally to the cell nucleus and e

xhibits both nuclear and plasma membrane chlorid  
e ion channel activity.