

小鼠抗 YBX2 单克隆抗体

中文名称：小鼠抗 YBX2 单克隆抗体

英文名称：Anti-YBX2 mouse monoclonal antibody

别名：Y-box binding protein 2; DBPC; MSY2; CSDA3; CONTRIN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：YBX2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes a nucleic acid binding protein which is highly expressed in germ cells. The encoded protein binds to a Y-box element in the promoters of certain genes but also binds to mRNA transcribed from these genes. Pseudo genes for this gene are located on chromosome 10 and 15.
Applications:	WB
Name of antibody:	YBX2
Immunogen:	Fusion protein of human YBX2
Full name:	Y-box binding protein 2
Synonyms:	DBPC; MSY2; CSDA3; CONTRIN
SwissProt:	Q9Y2T7
WB Predicted band size:	39 KD



WB Positive control:	Mouse testis tissue lysate
WB Recommended dilution:	500-2000