

小鼠抗 ZNF572 单克隆抗体

中文名称：小鼠抗 ZNF572 单克隆抗体

英文名称：Anti-ZNF572 mouse monoclonal antibody

别名：zinc finger protein 572

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：ZNF572

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB
Name of antibody:	ZNF572
Immunogen:	Fusion protein of human ZNF572
Full name:	zinc finger protein 572
SwissProt:	Q7Z3I7
WB Predicted band size:	61 KD
WB Positive control:	HepG2 cell lysate
WB Recommended dilution:	200-1000