

## Cy3 标记链霉亲和素

中文名称: Cy3 标记链霉亲和素

英文名称: Cy3-conjugated Streptavidin

相关类别: 抗原

储 存: 冷冻(-20℃)

## 概述

Cy3-conjugated Streptavidin

## 技术规格

Full name	Streptavidin
Applications	2-5 μg/ml for most applications
Background	Streptavidin is a tetrameric protein composed of identic subunits. The preparation contains an N- and C-terminal shortened variant (core streptavidin) with improved properties concerning homogen eity, solubility, resistance towards proteolytic degradation and acc essibility of the biotin binding pocket as compared to native stre ptavidin. The high affinity recognition of biotin and biotinylated molecules has made streptavidin one of the most important components in diagnostics and laboratory kits.