

RPS27L 抗原（重组蛋白）

中文名称： RPS27L 抗原（重组蛋白）

英文名称： RPS27L Antigen (Recombinant Protein)

别名： ribosomal protein S27 like

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 2-84 amino acids of human RPS27L

技术规格

Full name:	ribosomal protein S27 like
Swissprot:	Q71UM5
Gene Accession:	BC003667
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27, which suggests the encoded protein may be a component of the 40S ribosomal subunit.