

## RFLNA 抗原（重组蛋白）

中文名称： RFLNA 抗原（重组蛋白）

英文名称： RFLNA Antigen (Recombinant Protein)

别名： refilin A; CFM2; FAM101A

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 82-216 amino acids of human RFLNA

### 技术规格

<b>Full name:</b>	refilin A
<b>Synonyms:</b>	CFM2; FAM101A
<b>Swissprot:</b>	Q6ZTI6
<b>Gene Accession:</b>	BC141805
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Involved in the regulation of the perinuclear actin network and nuclear shape through interaction with filamins. Plays an essential role in actin cytoskeleton formation in developing cartilaginous cells.

**mlbio** 梅联生物  
Good elisakit producers