

## RGS6 抗原(重组蛋白)

中文名称: RGS6 抗原(重组蛋白)

英文名称: RGS6 Antigen (Recombinant Protein)

别 名: regulator of G protein signaling 6; GAP; S914; HA117

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to a region derived from 223-472 amino acids of human RGS6

## 技术规格

Full name:	regulator of G protein signaling 6
Synonyms:	GAP; S914; HA117
Swissprot:	P49758
Gene Accession:	BC113572
Purity:	>85%, as determined by Coomassie blue stained SDS-PA GE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the RGS (regulator of G protein signaling) family of proteins, which are defined by the presence of a RGS domain that confers the GTPase-a ctivating activity of these proteins toward certain G alpha subunits. This protein also belongs to a subfamily of RGS proteins characterized by the presence of DEP and GGL d



omains, the latter a G beta 5-interacting domain. The RG S proteins negatively regulate G protein signaling, and m ay modulate neuronal, cardiovascular, lymphocytic activitie s, and cancer risk. Many alternatively spliced transcript variants encoding different isoforms with long or short N-ter minal domains, complete or incomplete GGL domains, and distinct C-terminal domains, have been described for this gene, however, the full-length nature of some of these variants is not known.