

RGS11 抗原(重组蛋白)

中文名称: RGS11 抗原(重组蛋白)

英文名称: RGS11 Antigen (Recombinant Protein)

别 名: regulator of G-protein signaling 11; RS11

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to C terminal 250 amino acids of human RGS11

技术规格

Full name:	regulator of G-protein signaling 11				
Synonyms:	RS11				
Swissprot:	O94810				
Gene Accession:	BC141803				
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE				
Expression system:	Escherichia coli				
Tags:	His tag C-Terminus, GST tag N-Terminus				
Background:	The protein encoded by this gene belongs to the RGS (regul ator of G protein signaling) family. Members of the RGS family act as GTPase-activating proteins on the alpha subunits of heterotrimeric, signal-transducing G proteins. This protein inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits, thereby driving them into their inactive GDP-bound form. Alternative splicing occurs at this locu				



s and	four	transcript	variants	encoding	distinct	isoforms	have
been	identi	ified.					