

RMC1 抗原（重组蛋白）

中文名称： RMC1 抗原（重组蛋白）

英文名称： RMC1 Antigen (Recombinant Protein)

别名： regulator of MON1-CCZ1; MIC1; Mic-1; WDR98; C18orf8; HsT2591

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 458-657 amino acids of human RMC1

技术规格

Full name:	regulator of MON1-CCZ1
Synonyms:	MIC1; Mic-1; WDR98; C18orf8; HsT2591
Swissprot:	Q96DM3
Gene Accession:	BC002950
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a colon cancer associated protein. [provided by RefSeq, Jan 2013]