

RMND5B 抗原（重组蛋白）

中文名称： RMND5B 抗原（重组蛋白）

英文名称： RMND5B Antigen (Recombinant Protein)

别名： required for meiotic nuclear division 5 homolog B; GID2; GID2B

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 194-393 amino acids of human RMND5B

技术规格

Full name:	required for meiotic nuclear division 5 homolog B
Synonyms:	GID2; GID2B
Swissprot:	Q96G75
Gene Accession:	BC009911
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein remains unknown.