

## RNF128 抗原(重组蛋白)

中文名称: RNF128 抗原(重组蛋白)

英文名称: RNF128 Antigen (Recombinant Protein)

别 名: ring finger protein 128; GRAIL

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to a region derived from 39-207 amino acids of human RNF128

## 技术规格

Full name:	ring finger protein 128
Synonyms:	GRAIL
Swissprot:	Q8TEB7
Gene Accession:	BC063404
Purity:	>85%, as determined by Coomassie blue stained SDS-PAG E
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a type I transmembr ane protein that localizes to the endocytic pathway. This protein contains a RING zinc-finger motif and has been sh own to possess E3 ubiquitin ligase activity. Expression of t his gene in retrovirally transduced T cell hybridoma significantly inhibits activation-induced IL2 and IL4 cytokine pro



duction. Induced expression of this gene was observed in
anergic CD4(+) T cells, which suggested a role in the indu
ction of anergic phenotype. Alternatively spliced transcript
variants encoding distinct isoforms have been reported.