

## RNF185 抗原（重组蛋白）

中文名称： RNF185 抗原（重组蛋白）

英文名称： RNF185 Antigen (Recombinant Protein)

别名： ring finger protein 185

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to N terminal 185 amino acids of human RNF185

### 技术规格

<b>Full name:</b>	ring finger protein 185
<b>Swissprot:</b>	Q96GF1
<b>Gene Accession:</b>	BC009504
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	RNF185 (ring finger protein 185), also known as FLJ38628, is a 192 amino acid multi-pass membrane protein containing one RING-type zinc finger. Two RNF185 isoforms exist as a result of alternative splicing, and the gene encoding RNF185 maps to human chromosome 22, which houses over 500 genes and is the second smallest human chromosome. Mutations in several of the genes that map to chromosome 22 are involved in the development of Phelan-McDermid syndrome, Neurofibromatosis type 2, autism and schizophrenia.

**mlbio** 梅联生物  
Good elisakit producers