

RUFY4 抗原（重组蛋白）

中文名称： RUFY4 抗原（重组蛋白）

英文名称： RUFY4 Antigen (Recombinant Protein)

别名： RUN and FYVE domain containing 4; ZFYVE31

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 372-571 amino acids of human RUFY4

技术规格

Full name:	RUN and FYVE domain containing 4
Synonyms:	ZFYVE31
Swissprot:	Q6ZNE9
Gene Accession:	BC113712
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein remains unknown.