

## RWDD1 抗原（重组蛋白）

中文名称： RWDD1 抗原（重组蛋白）

英文名称： RWDD1 Antigen (Recombinant Protein)

别名： RWD domain containing 1; CGI-24; PTD013

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 1-243 amino acids of human RWDD1

### 技术规格

<b>Full name:</b>	RWD domain containing 1
<b>Synonyms:</b>	CGI-24; PTD013
<b>Swissprot:</b>	Q9H446
<b>Gene Accession:</b>	BC015802
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Protects DRG2 from proteolytic degradation.