

## SAMM50 抗原（重组蛋白）

中文名称： SAMM50 抗原（重组蛋白）

英文名称： SAMM50 Antigen (Recombinant Protein)

别名： SAMM50 sorting and assembly machinery component; OMP85; SAM50; TOB55; TRG-3; CGI-51; YNL026W

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 270-469 amino acids of human SAMM50

### 技术规格

<b>Full name:</b>	SAMM50 sorting and assembly machinery component
<b>Synonyms:</b>	OMP85; SAM50; TOB55; TRG-3; CGI-51; YNL026W
<b>Swissprot:</b>	Q9Y512
<b>Gene Accession:</b>	BC011681
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a component of the Sorting and Assembly Machinery (SAM) of the mitochondrial outer membrane. The Sam complex functions in the assembly of beta-barrel proteins into the outer mitochondrial membrane.