

SAP130 抗原（重组蛋白）

中文名称： SAP130 抗原（重组蛋白）

英文名称： SAP130 Antigen (Recombinant Protein)

别名： Sin3A associated protein 130

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 849-1048 amino acids of human SAP130

技术规格

Full name:	Sin3A associated protein 130
Swissprot:	Q9H0E3
Gene Accession:	BC017453
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	SAP130 is a subunit of the histone deacetylase (see HDAC1; MIM 601241)-dependent SIN3A (MIM 607776) corepressor complex (Fleischer et al., 2003 [PubMed 12724404]).