

## SASH1 抗原

- 中文名称: SASH1 抗原
- 英文名称: SASH1 Antigen
- 相关类别: 抗原
- 储 存: 冷冻 (-20℃)

## 概述

Synthetic peptide corresponding to a region derived from 78-92 amino acids of human SAM and SH3 domain containing 1

## 技术规格

Full name	SAM and SH3 domain containing 1
Synonyms	SH3D6A; dJ323M4; dJ323M4.1
Swissprot	O94885
Gene Accession	NP_056093
Background	SASH1 (SAM and SH3 domain-containing protein 1), also kn own as PEPE1 (Proline-glutamate repeat-containing protein), is a 1247 amino acid protein that is significantly downregul ated in the majority of primary breast tumor tissues, breast cancer cell lines, lung and thyroid tumors, as well as in cert ain colon carcinomas. It has been hypothesized that its expr ession is suppressed not due to mutation of the SASH1 ge ne, but instead via other mechanisms, such as promoter me thylation.

