

## SDHA 抗原(重组蛋白)

中文名称: SDHA 抗原(重组蛋白)

英文名称: SDHA Antigen (Recombinant Protein)

别 名: succinate dehydrogenase complex, subunit A, flavoprotein (Fp); FP; PGL5; SDH1;

SDH2; SDHF; CMD1GG

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to C terminal 215 amino acids of human SDHA

## 技术规格

Full name:	succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
Synonyms:	FP; PGL5; SDH1; SDH2; SDHF; CMD1GG
Swissprot:	P31040
Gene Accession:	BC001380
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a major catalytic subunit of succinate-ubiqui none oxidoreductase, a complex of the mitochondrial respirator y chain. The complex is composed of four nuclear-encoded su bunits and is localized in the mitochondrial inner membrane.  Mutations in this gene have been associated with a form of mi tochondrial respiratory chain deficiency known as Leigh Syndro



me. A pseudogene has been identified on chromosome 3q29.