

SELENOF 抗原(重组蛋白)

中文名称: SELENOF 抗原(重组蛋白)

英文名称: SELENOF Antigen (Recombinant Protein)

别 名: selenoprotein F; SEP15

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Fusion protein corresponding to a region derived from 97-165 amino acids of human SELENOF

技术规格

Full name:	selenoprotein F
Synonyms:	SEP15
Swissprot:	O60613
Gene Accession:	BC005294
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene belongs to the SEP15/seleno protein M family. The exact function of this protein is not kno wn; however, it has been found to associate with UDP-glucose: glycoprotein glucosyltransferase (UGTR), an endoplasmic reticul um(ER)-resident protein, which is involved in the quality control of protein folding. The association with UGTR retains this prote in in the ER, where it may play a role in protein folding. It has



also been suggested to have a role in cancer etiology. This pro tein is a selenoprotein, containing the rare amino acid selenocy steine (Sec). Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein m RNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop si gnal. Alternatively spliced transcript variants have been found f or this gene.