

SNAP47 抗原（重组蛋白）

中文名称： SNAP47 抗原（重组蛋白）

英文名称： SNAP47 Antigen (Recombinant Protein)

别名： synaptosome associated protein 47; SVAP1; HEL170; SNAP-47; C1orf142; ESFI5812; HEL-S-290

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-464 amino acids of human SNAP47

技术规格

Full name:	synaptosome associated protein 47
Synonyms:	SVAP1; HEL170; SNAP-47; C1orf142; ESFI5812; HEL-S-290
Swissprot:	Q5SQN1
Gene Accession:	BC018760
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Plays a role in intracellular membrane fusion.