

## SP20 细胞裂解液

中文名称: SP20 细胞裂解液

英文名称: SP20 Cell Lysate

别 名: SP20 细胞裂解液

相关类别: 抗原

储 存: 冷冻(-20℃) 避光

## 概述

SP20 cell lysates from mouse myeloma cells were cultured in the media containing 10% FBS.

## 技术规格

Backg round Cells were lysed in lysis buffer (50mM Tris(PH7.4), 150mM NaCl, 1% NP-40, 1% SD S). 100  $\mu$ g of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer.