

SRGAP3 抗原(重组蛋白)

- 中文名称: SRGAP3 抗原(重组蛋白)
- 英文名称: SRGAP3 Antigen (Recombinant Protein)
- 别名: SLIT-ROBO Rho GTPase activating protein 3; WRP; MEGAP; SRGAP2; ARHGAP14
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Fusion protein corresponding to N terminal 250 amino acids of human SRGAP3

技术规格

Full name:	SLIT-ROBO Rho GTPase activating protein 3
Synonyms:	WRP; MEGAP; SRGAP2; ARHGAP14
Swissprot:	O43295
Gene Accession:	BC039300
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	SRGAP3 is a 1099 amino acid protein containing 1 FCH domain, 1 Rho-GAP domain and 1 SH3 domain. Expressed highly in brai n, and in lower levels in kidney, WRP is thought to play a role in cell migration through its interaction with Cdc42 and Rac1. C dc42 and Rac1 are two intracellular signaling proteins that regul ate the multistep cell migration process. WRP downregulates Cd c42 and Rac1 activity, thereby impairing actin and microtubule



dynamics, the formation of protrusions, and total cell migration.
Defects in the gene encoding WRP have been linked to severe
idiopathic mental retardation. Three isoforms of WRP exist as a
result of alternative splicing events.