

STX7 抗原（重组蛋白）

中文名称： STX7 抗原（重组蛋白）

英文名称： STX7 Antigen (Recombinant Protein)

别名： syntaxin 7

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 131 amino acids of human STX7

技术规格

Full name:	syntaxin 7
Swissprot:	O15400
Gene Accession:	BC011975
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in protein trafficking from the plasma membrane to the early endosome (EE) as well as in homotypic fusion of endocytic organelles. Mediates the endocytic trafficking from early endosomes to late endosomes and lysosomes.