

TACSTD2 抗原（重组蛋白）

中文名称： TACSTD2 抗原（重组蛋白）

英文名称： TACSTD2 Antigen (Recombinant Protein)

别名： tumor-associated calcium signal transducer 2; EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 27-323 amino acids of human TACSTD2

技术规格

Full name:	tumor-associated calcium signal transducer 2
Synonyms:	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Swissprot:	P09758
Gene Accession:	BC009409
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.