

THADA 抗原（重组蛋白）

中文名称：THADA 抗原（重组蛋白）

英文名称： THADA Antigen (Recombinant Protein)

别名： THADA, armadillo repeat containing; GITA; ARMC13

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1754-1953 amino acids of human THADA

技术规格

Full name:	THADA, armadillo repeat containing
Synonyms:	GITA; ARMC13
Swissprot:	Q6YHU6
Gene Accession:	BC025773
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is the target of 2p21 chromosomal aberrations in benign thyroid adenomas. Single nucleotide polymorphisms (SNPs) in this gene may be associated with type 2 diabetes and polycystic ovary syndrome. The encoded protein is likely involved in the death receptor pathway and apoptosis.