

TKT 抗原（重组蛋白）

中文名称： TKT 抗原（重组蛋白）

英文名称： TKT Antigen (Recombinant Protein)

别名： transketolase; TK; TKT1; HEL107

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 200 amino acids of human TKT

技术规格

Full name:	transketolase
Synonyms:	TK; TKT1; HEL107
Swissprot:	P29401
Gene Accession:	BC009970
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a thiamine-dependent enzyme which plays a role in the channeling of excess sugar phosphates to glycolysis in the pentose phosphate pathway. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.