

TNFAIP3 抗原(重组蛋白)

中文名称: TNFAIP3 抗原(重组蛋白)

英文名称: TNFAIP3 Antigen (Recombinant Protein)

别 名: TNF alpha induced protein 3; A20; AISBL; OTUD7C; TNFA1P2

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 580-790 amino acids of human TNFAIP3

技术规格

Full name:	TNF alpha induced protein 3
Synonyms:	A20; AISBL; OTUD7C; TNFA1P2
Swissprot:	P21580
Gene Accession:	BC114480
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene was identified as a gene whose expression is rapidly induced by the tumor necrosis factor (TNF). The protein encod ed by this gene is a zinc finger protein and ubiqitin-editing en zyme, and has been shown to inhibit NF-kappa B activation as well as TNF-mediated apoptosis. The encoded protein, which h as both ubiquitin ligase and deubiquitinase activities, is involve d in the cytokine-mediated immune and inflammatory respons es. Several transcript variants encoding the same protein have



been found for this gene.