

TRAPPC4 抗原（重组蛋白）

中文名称： TRAPPC4 抗原（重组蛋白）

英文名称： TRAPPC4 Antigen (Recombinant Protein)

别名： trafficking protein particle complex 4; SBDN; TRS23; PTD009; CGI-104; HSPC172; SYNBINDIN

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 1-219 amino acids of human TRAPPC4

技术规格

Full name:	trafficking protein particle complex 4
Synonyms:	SBDN; TRS23; PTD009; CGI-104; HSPC172; SYNBINDIN
Swissprot:	Q9Y296
Gene Accession:	BC010866
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May play a role in vesicular transport from endoplasmic reticulum to Golgi.

mlbio 梅联生物
Good elisakit producers