

## TRIM26 抗原(重组蛋白)

中文名称: TRIM26 抗原(重组蛋白)

英文名称: TRIM26 Antigen (Recombinant Protein)

别 名: tripartite motif containing 26; AFP; RNF95; ZNF173

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to a region derived from 1-200 amino acids of human TRIM26

## 技术规格

Full name:	tripartite motif containing 26
Synonyms:	AFP; RNF95; ZNF173
Swissprot:	Q12899
Gene Accession:	BC024039
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the tripar tite motif (TRIM) family. The TRIM motif includes three zincbinding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplas mic bodies. Although the function of the protein is unknown, the RING domain suggests that the protein may have DNA -binding activity. The gene localizes to the major histocompa



tibility complex (MHC) class I region on chromosome 6. Alte
rnatively spliced transcript variants encoding the same protei
n have been found for this gene.