

TRMT1 抗原（重组蛋白）

中文名称：TRMT1 抗原（重组蛋白）

英文名称：TRMT1 Antigen (Recombinant Protein)

别名：tRNA methyltransferase 1; TRM1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 460-659 amino acids of human TRMT1

技术规格

Full name:	tRNA methyltransferase 1
Synonyms:	TRM1
Swissprot:	Q9NXH9
Gene Accession:	BC002492
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.