

TUBA1B 抗原（重组蛋白）

中文名称： TUBA1B 抗原（重组蛋白）

英文名称： TUBA1B Antigen (Recombinant Protein)

别名： tubulin alpha 1b; K-ALPHA-1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 252-451 amino acids of human TUBA1B

技术规格

Full name:	tubulin alpha 1b
Synonyms:	K-ALPHA-1
Swissprot:	P68363
Gene Accession:	BC017004
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.