

TUBA1C 抗原（重组蛋白）

中文名称： TUBA1C 抗原（重组蛋白）

英文名称： TUBA1C Antigen (Recombinant Protein)

别名： tubulin, alpha 1c; TUBA6; bcm948

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 200 amino acids of human TUBA1C

技术规格

Full name:	tubulin, alpha 1c
Synonyms:	TUBA6; bcm948
Swissprot:	Q9BQE3
Gene Accession:	BC004949
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.