

## TXLNA 抗原（重组蛋白）

中文名称： TXLNA 抗原（重组蛋白）

英文名称： TXLNA Antigen (Recombinant Protein)

别名： taxilin alpha; IL14; TXLN

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 186-412 amino acids of human TXLNA

### 技术规格

<b>Full name:</b>	taxilin alpha
<b>Synonyms:</b>	IL14; TXLN
<b>Swissprot:</b>	P40222
<b>Gene Accession:</b>	BC103823
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	$\alpha$ -taxilin is a novel binding partner of the syntaxin family which is implicated in intracellular vesicle trafficking. Through its C-terminal coiled-coil region, $\alpha$ -taxilin interacts with the nascent poly peptide-associated complex (NAC), which acts as a transcriptional coactivator.