

## TYK2 抗原(重组蛋白)

中文名称: TYK2 抗原(重组蛋白)

英文名称: TYK2 Antigen (Recombinant Protein)

别 名: tyrosine kinase 2; JTK1; IMD35

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to a region derived from 988-1187 amino acids of human TYK2

## 技术规格

Full name:	tyrosine kinase 2
Synonyms:	JTK1; IMD35
Swissprot:	P29597
Gene Accession:	BC014243
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the tyrosine kinase and, m ore specifically, the Janus kinases (JAKs) protein families. Thi s protein associates with the cytoplasmic domain of type I and type II cytokine receptors and promulgate cytokine sig nals by phosphorylating receptor subunits. It is also component of both the type I and type III interferon signaling pat hways. As such, it may play a role in anti-viral immunity. A



mutation in this gene has been associated with hyperimmu
noglobulin E syndrome (HIES) - a primary immunodeficiency
characterized by elevated serum immunoglobulin E.