

UBA5 抗原（重组蛋白）

中文名称： UBA5 抗原（重组蛋白）

英文名称： UBA5 Antigen (Recombinant Protein)

别名： ubiquitin-like modifier activating enzyme 5; THIFP1; UBE1DC1

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 205-404 amino acids of human UBA5

技术规格

Full name:	ubiquitin-like modifier activating enzyme 5
Synonyms:	THIFP1; UBE1DC1
Swissprot:	Q9GZZ9
Gene Accession:	BC009737
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the E1-like ubiquitin-activating enzyme family. This protein activates ubiquitin-fold modifier 1, a ubiquitin-like post-translational modifier protein, via the formation of a high-energy thioester bond. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. A pseudogene located on chromosome 1 has also been identified.

mlbio 梅联生物
Good elisakit producers