

## UBE2A 抗原（重组蛋白）

中文名称： UBE2A 抗原（重组蛋白）

英文名称： UBE2A Antigen (Recombinant Protein)

别名： ubiquitin conjugating enzyme E2 A; UBC2; HHR6A; MRXSN; RAD6A; MRXS30

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 20-152 amino acids of human UBE2A

### 技术规格

<b>Full name:</b>	ubiquitin conjugating enzyme E2 A
<b>Synonyms:</b>	UBC2; HHR6A; MRXSN; RAD6A; MRXS30
<b>Swissprot:</b>	P49459
<b>Gene Accession:</b>	BC010175
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, ubiquitin-conjugating enzymes, and ubiquitin-protein ligases. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme

is required for post-replicative DNA damage repair, and may play a role in transcriptional regulation. Mutations in this gene are associated with mental retardation. Alternative splicing results in multiple transcript variants.